



Humanware™





Explorē Product Line



Explorē 3





Explorē 5




Explorē 8

Explorē 8 Overview




Setup the Device

- Press&Hold the contrast button(Orange button)
- Setup the users preferred zoom
- Setup the enabled contrast
- Setup the text size
- Setup the language




Live Reading Mode

- Zoom in /out
- Change contrast
- Freeze picture
- Save picture



Distance Mode

- Zoom in /out
- Change contrast
- Freeze picture
- Save picture



Gallery

- View saved images
 - Using swipe gesture
 - Using arrows
- Change contrast
- Select an image
- Zoom in/Out
- Pan



Explorē 8 Key Features

- **Image Quality** → True and vivid colors with edge enhancements
- **Big Screen** (8 inches) → Bigger viewing area, a third more compared to the 7 inch.
- **Touchscreen** → Simplicity of use. Same gestures as on a tablet or smartphone.
- **Quick Startup Time** → Ready to use at all times with a click of a button.
- **Two cameras** → One 21Mpx camera for distance and one 21Mpx camera for close-up.
- **Contrast Colors** → 16 contrast colors and selection is customizable to user's preference.
- **4,5 hours Battery**

explorē 8 Key Features

- **Light weight (580g) and Small size** → Bring it with you anywhere.
- **Quick USB Charge** → USB C charging capabilities, 2 hours
- **Intuitive** → Control buttons below the screen for ease of use
- **Polarized LEDs** → Less reflection from glossy paper like magazines
- **Sturdy reading stand** → When folded (Distance) or unfolded (Close-up) triggers a camera switch automatically.
- **Accessories** → c-type USB cable, HDMI cable, protective membrane, cleaning cloth and polyurethane carrying case.

explorē 8 - What's New in 1.1.0

- Reduced drag effect when moving the device on a text.
- New auto focus for the distance camera by offering Manual and Auto modes in the Menu.
- Less refocus triggering in distance viewing auto mode.
- Touch to focus feature in distance viewing manual mode.
- Lights on/off control available for distance camera.